

QH-1


Technology: Axis Knee System

Description: Axis[™] Knee System is an affordable, high-quality and FDA-approved total knee replacement system specifically designed for the Asian population. Its instrumentation and surgical technique allow surgeons to perform the procedure accurately, eliminating the possibility of implant failure.

Contact Details:

Engr. Jude Sasing
Orthopaedic International Inc.
0917 887 0256





QH-2

Technology: eHATID

Description: eHealth Tablet for Informed Decision Making of LGUs (eHATID) provides health information system (HIS) and decision-making support to LGUs through an Electronic Medical Record that generates reports for DOH and PhilHealth. It aims to improve efficiency and transparency using an enhanced HIS model that promotes health information convergence and informed decision-making.

Contact Details:

Dr. Dennis Batangan Ateneo de Manila University (02) 426-6067 locals 213 or 209



QH-3

Technology: RxBox

Description: RxBox is a telemedicine device for capturing medical signals via built-in medical sensors, storing data in an electronic medical record and transmitting health information via internet to a clinical specialist for expert advice. RxBox can reduce overall cost of healthcare by enabling health workers to diagnose, monitor and treat patients within the rural health facility.

Contact Details:

Dr. Portia F. Marcelo National Telehealth Center 0927 881 7929



QH-4

Technology: Ventilator

Description: It is a compact, easy-to-use, volume-cycled servo ventilator. It is AC/DC-powered with backup battery, and it can also be used as transport and home ventilator. Its embedded software incorporates a self-diagnosis protocol for offsite malfunction troubleshooting.

Contact Details:

Dr. Abundio Balgos

University of the Philippines Manila

0917 848 5642



QH-5

Technology: OL Trap

Description: Ovicidal/Larvicidal Trap (OL Trap) is used to control the population of dengue-carrying Aedes aegypti mosquito. It works by attracting female mosquitoes to lay eggs on the paddle drenched with an organic solution which kills the eggs and larvae, preventing mosquitoes from reaching adulthood. It is potent for mosquitoes but safe for humans

Contact Details:

Dr. Maria Patricia V. Azanza Industrial Technology Development Institute (02) 837 3167



QH-6

Technology: Biotek-M

Description: Biotek-M kit is used as a highly efficient confirmatory test for diagnosis of dengue infection within first 5 days of illness, with results in an hour or less. Using Biotek-M would mean less admission for denguesuspected cases, therefore saving resources for both hospital and patients.

Contact Details:

Dr. Raul V. Destura University of the Philippines Manila (02) 526 4349

